


Applications (CUY21)

● Electroporation-mediated gene transfer system applied to cultured CNS neurons [Publication 14]




Schematic representation of an electroporation set-up.

A fragment of the mouse embryonic hippocampus was placed on a Millipore membrane filter and 5 μ l EP buffer containing 1mg/ml of plasmid DNA was applied onto the tissue.

A tungsten needle was attached to the surface of a droplet.

After application of square pulses the tissue fragment was returned to a petri dish containing ice-cold HBSS solution.



Electroporation-mediated expression of fluorescent proteins in hippocampal neurons.

(a-c) Organ culture of hippocampal tissue fragments three days after electroporation with CAG-eGFP (a), Ta1X4-eGFP (b), and β -actin-eGFP (c) expression constructs.

(d, e) A mature hippocampal neuron maintained 14 days in dissociated culture after electroporation of a β -actin-eGFP expression construct. Higher magnification view of the region marked by a rectangle in (d) reveals dendritic spines on the surface of dendritic shafts (arrows in e).

(f, g) A hippocampal neuron 7 days after electroporation of 1:1 mixture of Ta1X4-eGFP and Ta1X4-mRFP1. Both eGFP fluorescence (f) and mRFP1 fluorescence (g) can be observed in a single cell.

(h) Relative fluorescence intensity of hippocampal tissue fragments after electroporation of eGFP-expression plasmids with four different promoter sequences. The tissue fragments were maintained in culture for 4 days, fixed and observed under a confocal microscope. Fluorescence intensities per unit area of the tissue fragments were determined.

(i) Relative fluorescence intensity of hippocampal tissue fragments isolated at two different developmental stages and electroporated with β -actin-eGFP. Tissue fragments were maintained for 4 days in culture and subsequently fixed. Fluorescence intensities were measured using a confocal microscope.

Bars = 50 μ m (a-d, f, g); 10 μ m (e).