

# NEPA GENE Cell Transfection Data

## Primary Cells: With Electroporation Cuvettes

V:Viability / TE:Transfection Efficiency

Cell Line	Description, Species	Characteristics, etc.	V	TE	IMG
HUVEC	<i>Human Umbilical Vein Endothelial cells</i>	Primary, Normal, Adherent	95%	75%	
HASM	<i>Human Airway Smooth Muscle cells</i>	Primary, Adherent	90%	80%	Y
	<i>Human Endometrial Stromal cells</i>	Primary Adherent	95%	90%	Y
	<i>Human Uterine Cervical Fibroblasts</i>	Primary, Adherent	65%	90%	Y
	<i>Human Dermal Fibroblasts</i>	Primary, Adherent	95%	89%	Y
	<i>Human Malignant Mesothelioma cells</i>	Primary, Adherent	75%	55%	
	<i>Human T cells</i>	Primary, Suspension	58%	90%	Y
	<i>Human T cells</i>	Primary, Suspension	50%	76%	
	<i>Human NK cells</i>	Primary, Suspension	31%	47%	Y
PBMC	<i>Peripheral Blood Mononuclear cells</i>	Primary, Suspension	93%	66%	
	<i>Human Chronic Lymphocytic Leukemia (CLL)</i>	Primary, Suspension	82%	70%	Y
	<i>Mouse Cerebral Cortex Neurons (Embryonic day14)</i>	Primary, Adherent	80%	70%	Y
	<i>Mouse Hippocampal Neurons (Embryonic day 17)</i>	Primary, Adherent	65%	70%	Y
	<i>Mouse Hippocampal Neurons (Embryonic day 14)</i>	Primary, Adherent	80%	60%	
	<i>Mouse Neural Progenitor Cells</i>	Primary, Adherent	80%	60%	Y
	<i>Mouse Basal Ganglia Primordium (Embryonic day 13.5)</i>	Primary, Adherent	91%	71%	

	<i>Mouse Cerebellar Granule Neurons</i>	Primary, Adherent	91%	65%	Y
	<i>Mouse DRG Neurons</i>	Primary, Adherent	70%	70%	Y
MEF	<i>Mouse Embryonic Fibroblasts</i>	Primary, Adherent	90%	85%	Y
MEF	<i>Mouse Embryonic Fibroblasts</i>	Primary, Adherent	75%	85%	Y
	<i>Mouse Embryonic Skin Fibroblasts</i>	Primary, Adherent	80%	50%	Y
	<i>Mouse Cervical Epithelial cells</i>	Primary, Adherent	82%	55%	
	<i>Mouse Vascular Adventitial Fibroblasts</i>	Primary, Adherent	90%	50%	Y
BMMC (BMMCs)	<i>Mouse Bone Marrow-Derived Mast cells</i>	Primary, Suspension	64%	78%	Y
BMMC (BMMCs)	<i>Mouse Bone Marrow-Derived Mast cells</i>	Primary, Suspension	80%	83%	Y
	<i>Mouse Liver cells</i>	Primary, Adherent	75%	65%	Y
	<i>Mouse Osteoblast cells</i>	Primary, Adherent	85%	60%	
	<i>Mouse Muscle cells (extensor digitorum longus)</i>	Primary, Adherent	68%	54%	Y
	<i>Rat Cerebral Cortex Neurons (Embryonic day16)</i>	Primary, Adherent	70%	75%	Y
	<i>Rat Hippocampal Neurons</i>	Primary, Adherent	60%	80%	
	<i>Rat Bulbar Neurons (Embryonic day 15)</i>	Primary, Adherent	80%	75%	Y
	<i>Rat Cerebellar Neurons</i>	Primary, Adherent	70%	55%	
	<i>Rat Cerebellar Granule cells</i>	Primary, Adherent	70%	80%	
	<i>Rat Schwann cells</i>	Primary, Adherent	90%	80%	Y
	<i>Rat Schwann cells</i>	Primary, Adherent	90%	60%	Y
	<i>Rat Müller cells</i>	Primary, Adherent	90%	50%	Y

	<i>Rat Meningeal Fibroblasts</i> (Postnatal day 3)	Primary, Adherent	90%	95%	Y
OEC	<i>Rat Olfactory Ensheathing cells</i> (Postnatal week 3)	Primary, Adherent	93%	46%	Y
REF	<i>Rat Embryonic Fibroblasts</i>	Primary, Adherent	65%	65%	Y
	<i>Chick Embryonic Fibroblasts</i>	Primary Adherent	80%	90%	Y
	<i>Chick Embryonic Cerebellar Granule Cells</i> (Embryonic day 11)	Primary, Adherent	86%	83%	
	<i>Rabbit Spleen cells (B cells)</i>	Primary, Suspension	70%	45%	Y
	<i>Goat Embryonic Epithelial Fibroblasts</i> (Embryonic day 43)	Primary, Adherent	80%	55%	Y

## Stem Cells, ES Cells and iPS Cells: With Electroporation Cuvettes

Cell Line	Description, Species	Characteristics, etc.	V	TE	IMG
	<i>Human iPS cells</i>	Adherent	Excellent	Excellent	Y
	<i>Human iPS cells (201B7)</i>	Adherent	80%	69%	
	<i>Human iPS cells</i> After EP: Feeder Free Culture	Adherent		73%	Y
	<i>Human ES cells</i>	Adherent	Excellent	Excellent	Y
	<i>Human ES cells (H9 p.51)</i>	Adherent	50-60%	50-60%	Y
	<i>Human Mesenchymal Stem cells</i>	Primary, Cord blood derived	78%	75%	Y
	<i>Human Mesenchymal Stem cells</i>	Bone Marrow derived	70%	80%	Y
TS	<i>Mouse Trophoblast Stem cells</i>	Adherent	59%	47%	Y
C3H/10T1/2	<i>Mouse Embryonic Mesenchymal Stem cells</i>	Adherent	70%	85%	Y
TT2	<i>Mouse TT2 ES cells</i>	Adherent	55%	55%	Y
	<i>Mouse ES cells, 129 strain, R1/E</i>	Adherent	80%	75%	Y

	<i>Mouse ES cells</i>	Adherent	80%	68%	Y
	<i>Mouse ES cells</i>	Adherent	74%	88%	Y
	<i>Mouse ES cells</i>	Adherent	80%	90%	Y
	<i>Mouse ES cells</i>	Adherent	70%	100%	Y
	<i>Mouse ES cells</i>	Adherent	77%	73%	
	<i>Mouse iPS cell derived Neural Stem cells</i>			86%	Y
	<i>Mouse Neural Stem cells</i>		90%	80%	Y
	<i>Mouse Neural Stem cells</i>		80%	60%	
	<i>Mouse Neurospheres</i> <i>The cells were from E13.5 mouse brain: ganglionic eminence</i>		90%	75%	Y
	<i>Mouse Neurospheres</i> <i>Neural Stem cells Derived from Mouse SVZ</i>		Excellent	Excellent	Y
	<i>Rat WDA ES-like cells</i>	Adherent	60%	80%	

## Cell Lines: With Electroporation Cuvettes

V: Viability / TE: Transfection Efficiency

Cell Line	Description, Species	Characteristics, etc.	V	TE	IMG
HeLa	<i>Human Cervical Carcinoma cells</i>	Epithelial, Immortalized Adherent	87%	93%	Y
HeLa-K	<i>Human Cervical Carcinoma cells</i>	Epithelial, Immortalized	90%	90%	
293 (HEK293)	<i>Human Embryonic Kidney cells</i>	Epithelial, Adenovirus type 5	92%	91%	Y
293 (HEK293)	<i>Human Embryonic Kidney cells</i>	Epithelial, Adenovirus type 5	90%	90%	
293 (HEK293)	<i>Human Embryonic Kidney cells</i>	Epithelial, Adherent	90%	70%	Y
293 (HEK293)	<i>Human Embryonic Kidney cells</i>	Epithelial, Adherent	72%	85%	Y

293T (HEK293T)	<i>Human Embryonic Kidney cells</i>	Epithelial, SV40 large T antigen	90%	95%	Y
293T (HEK293T)	<i>Human Embryonic Kidney cells</i>	Epithelial, SV40 large T antigen	83%	87%	Y
293T (HEK293T)	<i>Human Embryonic Kidney cells</i>	Epithelial, SV40 large T antigen	70%	99%	
TIG-3	<i>Human Embryonic Lung Fibroblasts</i>	Diploid, Normal Adherent	90%	80%	Y
TIG-7	<i>Human Embryonic Lung Fibroblasts</i>	Diploid, Normal Adherent	89%	76%	Y
MRC-5	<i>Human Embryonic Lung Fibroblasts</i>	Diploid, Normal Adherent	85%	90%	Y
WI-38	<i>Human Embryonic Lung Fibroblasts</i>	Adherent	80%	70%	
HDF	<i>Human Dermal Fibroblasts (106-05)</i>	Adherent	90%	90%	Y
HaCaT	<i>Human Keratinocyte cells</i>	Normal, Immortalized	96%	69%	Y
BEAS-2B	<i>Human Bronchial Epithelial cells</i>	Adherent	75%	96%	
SUSM-1	<i>Human Fibroblasts</i>	Fetal Liver, Epithelial,	77%	71%	Y
KMST-6	<i>Human Fibroblasts</i>	Epithelial Immortalized	70%	60%	Y
HT1080	<i>Human Fibrosarcoma cells</i>	Connective Tissue Acetabulum,	93%	81%	Y
HT1080	<i>Human Fibrosarcoma cells</i>	Connective Tissue Acetabulum,	80%	90%	
MG-63	<i>Human Osteosarcoma Cells</i>	Adherent	70%	80%	
U2OS	<i>Human Osteosarcoma Cells</i>	Adherent	70%	80%	
Saos-2	<i>Human Osteosarcoma Cells</i>	Adherent	60%	75%	
PANC-1	<i>Human Pancreatic Carcinoma cells</i>	Adherent	55%	75%	
MIA-PaCa-2	<i>Human Pancreatic Carcinoma cells</i>	Epithelial, Adherent	80%	77%	Y
HepG2	<i>Human Hepatoma cells</i>	Epithelial, Adherent	95%	85%	Y

HepG2	<i>Human Hepatoma cells</i>	Epithelial, Adherent	80%	76%	Y
HuH-7	<i>Human Hepatoma cells</i>	Epithelial, Adherent	82%	85%	
TFK-1	<i>Human Bile Duct Adenocarcinoma cells</i>	Adherent	50%	70%	Y
H69	<i>Human Small-Cell Lung Cancer cells</i>	Adhereht	90%	85%	Y
H1299 (NCI-H1299)	<i>Human Lung Cancer cells</i>	Non-Small, Epithelial, Adherent	90%	90%	
H1299	<i>Human Lung Cancer cells</i>	Non-Small, Epithelial	80%	80%	Y
KB31	<i>Human Epidermoid Carcinoma cells</i>	Adherent	58%	64%	Y
HSC-2	<i>Human Squamous Carcinoma cells</i>	Oral , Adherent	90-95%	98%	Y
HSC-3	<i>Human Squamous Carcinoma cells</i>	Tongue , Adherent	90-95%	98%	Y
Ca9-22	<i>Human Squamous Carcinoma cells</i>	Gingival, Adherent	60%	60%	Y
HGF	<i>Human Gingival Fibroblasts</i>	Normal, Adherent	Excellent	Excellent	Y
HEp-2	<i>Human Laryngeal Carcinoma cells</i>	Epithelial, Adherent	70%	90%	Y
MCF-7	<i>Human Breast Cancer cells</i>	Epithelial, Adherent	95%	80%	Y
MCF-7	<i>Human Breast Cancer cells</i>	Epithelial, Adherent	80%	70%	
MCF-7	<i>Human Breast Cancer cells</i>	Epithelial, Adherent	81%	65%	Y
T47D	<i>Human Breast Cancer cells</i>	Ductal Carcinoma Epithelial	90%	80-90%	
BT-20	<i>Human Breast Cancer cells</i>	Adherent	70%	80%	Y
MDA-MB-231	<i>Human Breast Cancer cells</i>	Adherent	85%	>90%	Y
MCF 10A	<i>Human Breast Cells</i>	Normal, Mammary Epithelial	90%	80%	Y
MCF 10A	<i>Human Breast Cells</i>	Normal, Mammary Epithelial	Excellent	Excellent	Y

NUGC-3	<i>Human Gastric Carcinoma cells</i>	Adherent	73%	68%	
A549	<i>Human Lung Adenocarcinoma cells</i>	Epithelial, Adherent	80-90%	90%	
TE-1	<i>Human Esophageal Carcinoma cells</i>	Adherent	80-90%	41%	Y
LNCaP	<i>Human Prostate Carcinoma</i>	Epithelial, Adherent	71%	90%	
PC-3	<i>Human Prostate Cancer cells</i>	Epithelial, Adherent	90%	95%	
PC-3	<i>Human Prostate Cancer cells</i>	Epithelial, Adherent	86%	55%	Y
PNT2	<i>Human Prostate Epithelial cells</i>	Epithelial, Adherent	85%	80%	Y
LoVo	<i>Human Colon Adenocarcinoma cells</i>	Epithelial, Adherent	85%	60%	Y
HCT116	<i>Human Colon Cancer cells</i>	Adherent	80%	95%	Y
HCT116	<i>Human Colon Cancer cells</i>	Adherent	95%	90%	Y
HCT116	<i>Human Colon Cancer cells</i>	Adherent	80%	80%	Y
Caco-2	<i>Human Colon Cancer cells</i>	Adherent	95%	80%	Y
SKOV-3	<i>Human Ovarian Carcinoma cells</i>	Adherent	90%	90%	Y
OVCAR-3	<i>Human Ovarian Carcinoma cells</i>	Adherent	90%	79%	Y
RMG-1	<i>Human Ovarian Clear Cell Adenocarcinoma</i>	Adherent	97%	67%	Y
SK-N-SH	<i>Human Neuroblastoma cells</i>	Epithelial, Adherent	95%	95%	Y
SH-SY5Y	<i>Human Neuroblastoma cells</i>	Adherent	60%	90%	Y
SH-SY5Y	<i>Human Neuroblastoma cells</i>	Adherent	70%	70%	
SH-SY5Y	<i>Human Neuroblastoma cells</i>	Adherent			Y
NB9	<i>Human Neuroblastoma cells</i>	Adherent	70%	70%	

NB69	<i>Human Neuroblastoma cells</i>	Adherent	95%	80%	
NB-39-nu	<i>Human Neuroblastoma cells</i>	Adherent	60%	63%	
KG-1-C	<i>Human Oligodendroglial cells</i>	Adherent	85%	60%	Y
A172	<i>Human Glioblastoma cells</i>	Adherent	80-90%	70%	Y
NP3	<i>Human Glioblastoma cells</i>	Adherent	98%	62%	Y
1321N1	<i>Human Astrocytoma cells</i>	Adherent	80%	80%	Y
	<i>Immortalized Human Pericytes</i>	Adherent	83%	50%	Y
iHAM-4	<i>Human Amniotic Mesenchymal cells</i>	Adherent	59%	95%	
	<i>Human Dental Pulp cells</i>	Adherent	90%	85%	Y
	<i>Human Dental Pulp cells</i>	Adherent	85%	69%	Y
HTR-8/Svneo	<i>Human Trophoblast cells</i>	Adherent	95%	67%	Y
SRA 01/04	<i>Human Lens Epithelial cells</i>	Adherent	97%	80%	Y
RPE	<i>Human Retinal Pigment Epithelium cells</i>	Adherent	90%	70%	
RPE-1 (hTERT RPE-1)	<i>Human Retinal Pigment Epithelium cells</i>	Adherent Immortalized	95%	80%	Y
RPTEC	<i>Human Renal Proximal Tubule Epithelial cells</i>	Adherent	70%	85%	
UT-7/TPO	<i>Human Megakaryocytic cells</i>	Suspension	51%	46%	Y
P493-6	<i>Human B cells</i>	Suspension	80%	40-50%	Y
MV4-11	<i>Human Acute Myeloid Leukemia cells</i>	Suspension	70%	60%	Y
Jurkat	<i>Human T-cell Leukemia cells</i>	Round Single, Suspension	99%	92%	Y
Jurkat	<i>Human T-cell Leukemia cells</i>	Round Single, Suspension	73%	94%	Y



Jurkat	<i>Human T-cell Leukemia cells</i>	Round Single, Suspension	85%	85%	Y
Jurkat	<i>Human T-cell Leukemia cells</i>	Round Single, Suspension	89%	85%	
Jurkat	<i>Human T-cell Leukemia cells</i>	Round Single, Suspension	95%	95%	
Hut78	<i>Human T-cell Lymphoma cells</i>	Suspension	51%	74%	
Jeko-1	<i>Human Mantle Cell Lymphoma (MCL) cells</i>	Peripheral blood, Suspension	80%	63%	
MOLT-4	<i>Human Acute Lymphoblastic Leukemia cells</i>	Suspension	95%	70%	
697	<i>Human Pre-B Acute Lymphoblastic Leukemia cells</i>	Suspension	68%	93%	
Nalm-6	<i>Human B-cell Precursor Leukemia cells</i>	Blood, Suspension	77%	82%	
Nalm-6	<i>Human B-cell Precursor Leukemia cells</i>	Blood, Suspension	67%	70%	
KG-1	<i>Human Acute Myeloid Leukemia cells</i>	Suspension	60%	65%	Y
KG-1	<i>Human Acute Myeloid Leukemia cells</i>	Suspension	70%	65%	
Kasumi-1	<i>Human Acute Myeloid Leukemia cells</i>	Suspension	66%	79%	
USCD/AML1	<i>Human Leukemia cells</i>	Suspension	50%	50%	
M7	<i>Human Acute Non Lymphocytic Leukemia cells</i>	Suspension	85%	80%	Y
KOPT-K1	<i>Human T cell acute lymphoblastic leukemia (T-ALL) cells</i>	Suspension	80%	60%	Y
	<i>Human EBV-immortalized B cells</i>	Suspension	58%	53%	Y
Namalwa	<i>Human Burkitt's Lymphoma cells</i>	B lymphocyte, Suspension	70%	75%	
Raji	<i>Human Burkitt's Lymphoma cells</i>	B lymphocyte, Suspension	97%	83%	Y
Toledo	<i>Human Burkitt's Lymphoma cells</i>	B lymphocyte, Suspension	81%	60%	Y
Mutu I	<i>Human Burkitt Lymphoma cells</i>	EBV Positive, Suspension	87%	91%	Y

Mutu III	<i>Human Burkitt Lymphoma cells</i>	EBV Positive, Suspension	54%	92%	
LCL	<i>Human Burkitt's Lymphoma cells</i>	Immortalized, Suspension	55%%	43%	Y
Ramos	<i>Human Burkitt's Lymphoma cells</i>	B lymphocyte, Suspension	83%	57%	
Ramos-Blue™	<i>Human Burkitt's Lymphoma cells</i>	B lymphocyte, Suspension	80%	55%	Y
BJAB	<i>Human EBV-negative Burkitt Lymphoma cells</i>	Suspension	96%	96%	Y
HL-60	<i>Human Promyelocytic Leukemia cells</i>	Suspension	80%	80%	
KU812	<i>Human Basophilic Leukaemia cells</i>	Suspension	96%	42%	
K562	<i>Human Chronic Myelogenous Leukemia cells</i>	Suspension	91%	99%	
K562	<i>Human Chronic Myelogenous Leukemiacells</i>	Suspension	77%	64%	
THP-1	<i>Human Acute Monocytic Leukemia cells</i>	Suspension	76%	63%	Y
THP-1	<i>Human Acute Monocytic Leukemia cells</i>	Suspension	56%	64%	Y
THP-1	<i>Human Acute Monocytic Leukemia cells</i>	Suspension	76%	59%	Y
THP-1	<i>Human Acute Monocytic Leukemia cells</i>	Suspension	85%	67%	
THP-1	<i>Human Acute Monocytic Leukemia cells</i>	Suspension	67%	85%	
KHYG-1	<i>Human Natural Killer (NK)Leukemia cells</i>	Suspension	51%	71%	Y

NIH/3T3	<i>Mouse Embryonic Fibroblasts</i>	Swiss NIH , Embryo	100%	90%	Y
NIH/3T3	<i>Mouse Embryonic Fibroblasts</i>	Swiss NIH , Embryo	74%	81%	Y
NIH/3T3	<i>Mouse Embryonic Fibroblasts</i>	Swiss NIH , Embryo	70%	70%	Y
PT67	<i>Mouse Fibroblasts (RetroPack PT67 cell line)</i>	Adherent	91%	66%	Y
3T3-L1	<i>Mouse Embryonic Fibroblasts (preadipocytes)</i>	Fibroblast like	90%	90%	Y

MEF	<i>Mouse Embryonic Fibroblasts</i>	Adherent	80%	90%	
MEF	<i>Mouse Embryonic Fibroblasts</i>	Adherent	75%	70%	Y
STO	<i>Mouse Embryonic fibroblasts</i>	Adherent	60%	51%	Y
N7	<i>Mouse Embryonic Hypothalamic cells (immortalized)</i>	Adherent	75%	100%	Y
P19C6	<i>Mouse Embryonic Carcinoma cells</i>	Adherent	90%	50%	Y
HL-1	<i>Mouse Cardiac Muscle cells</i>	Adherent	70%	70%	
L	<i>Mouse Fibroblasts</i>	Adherent	90%	60-70%	Y
B16	<i>Mouse Melanoma cells</i>	Adherent	70%	50%	Y
MC3T3-E1	<i>Mouse Osteoblastic cells</i>	Fibroblastic, Adherent	85%	75%	Y
C2C12	<i>Mouse Myoblast cells</i>	Adherent	94%	90%	Y
C2C12	<i>Mouse Myoblast cells</i>	Adherent	90%	90%	Y
C2C12	<i>Mouse Myoblast cells</i>	Adherent	80%	70%	Y
C2C12	<i>Mouse Myoblast cells</i>	Adherent	90%	95%	
NMuMG	<i>Mouse Mammary Gland Epithelial cells</i>	Adherent	80%	65%	
LLc1 (LL/2)	<i>Lewis Lung Cell Carcinoma 1 cells</i>	Adherent	87%	81%	Y
4T1	<i>Mouse Breast Cancer cells</i>	Adherent	90%	95%	Y
Colon-26	<i>Mouse Colon Adenocarcinoma cells</i>	Adherent	95%	90%	Y
Hepa1-6	<i>Mouse Hepatoma cells</i>	Adherent	50%	98%	
LM8	<i>Mouse Osteosarcoma cells</i>	Adherent	90%	85%	Y
AtT-20	<i>Mouse Pituitary Tumor cells</i>	Adherent	80%	80%	Y

MS-1	<i>Mouse Pancreatic Endothelial cells</i>	Adherent	90%	90%	
	<i>ddy Mouse Endometrial cells</i>	Fibroblast, Adherent	60%	80%	Y
Neuro-2a	<i>Mouse Neuroblastoma cells</i>	Adherent	90%	90%	Y
BV-2	<i>Mouse Microglial cells</i>	Adherent	65%	70%	Y
BV-2	<i>Mouse Microglial cells</i>	Adherent	90%	50%	Y
WR19L	<i>Mouse T-cell lymphoma cells</i>	Suspension	92%	60%	Y
BA/F3	<i>Mouse pro-B cells</i>	Suspension	90%	90%	
A20	<i>Mouse B-cell Lymphoma cells</i>	Suspension	70%	65%	Y
P815	<i>Mouse Mastocytoma cells</i>	Suspension	67%	68%	
J774.1	<i>Mouse Macrophage-like cells</i>	Monocyte-like	100%	70%	Y
RAW264.7	<i>Mouse Macrophage-like cells</i>	Monocyte-like, Adherent	70%	56%	Y
RAW264.7	<i>Mouse Macrophage-like cells</i>	Monocyte-like, Adherent	Excellent	Excellent	Y
RAW-D	<i>Mouse Macrophage cells</i>	Monocyte-like, Adherent	80%	80%	Y
MIN6	<i>Mouse Pancreatic Beta cells</i>	Suspension	57%	71%	Y
MEL	<i>Mouse Erythroleukemia cells</i>	Suspension	70%	50%	
XS106	<i>Mouse Dendritic Cells</i>	Suspension	61%	45%	Y
mDC	<i>Mouse Myeloid Dendritic cells</i>	Suspension	79%	72%	
416B	<i>Mouse Primitive Myeloid cells</i>	Suspension	89%	64%	
32D	<i>Mouse Myeloid cells</i>	Suspension		88%	
MC/9	<i>Mouse Mast cells</i>	Suspension	76%	84	Y

BMBa	<i>Mouse Bone marrow-derived basophils</i>	Suspension	45%	67%	Y
------	--	------------	-----	-----	---

PC12	<i>Rat Adrenal Pheochromocytoma cells</i>	Adherent	90%	70%	Y
H9c2	<i>Rat Ventricular Myoblasts</i>	Skeletal Muscle-like	71%	82%	Y
H9c2	<i>Rat Ventricular Myoblasts</i>	Skeletal Muscle-like	75%	80%	Y
REF	<i>Rat Embryonic Fibroblasts</i>	Adherent	90%	99%	Y
C6	<i>Rat Glioma cells</i>	Fibroblast-like, GFAP Positive,	80%	67%	Y
RSC96	<i>Rat Schwann cells</i>	Adherent	70%	85%	Y
TtT/GF	<i>Rat Anterior Pituitary cells</i>	Folliculostellate, Adherent	65%	83%	Y
UMR106	<i>Rat Osteoblastic cells</i>	Adherent	80%	70%	

CHO	<i>Chinese Hamster Ovary cells</i>	Epithelial, Adherent	74%	90%	Y
CHO-K1	<i>Chinese Hamster Ovary cells</i>	Epithelial, Adherent	95%	95%	Y
CHO-K1	<i>Chinese Hamster Ovary cells</i>	Epithelial, Adherent	90%	99%	
CHO-S	<i>Chinese Hamster Ovary cells</i>	Epithelial, Suspension	Excellent	56%	
COS-7	<i>African Green Monkey SV40-transf'd Kidney fibroblasts</i>	Adherent	61%	89%	
MDCK	<i>Madrin-Darby Canine Kidney cells</i>	Epithelial-like, Adherent	90%	95%	Y
BFF	<i>Bovine Fetal Fibroblasts</i>	Adherent	93%	71%	
BFF	<i>Bovine Fetal Fibroblasts</i>	Adherent	78%	72%	Y
	<i>Bovine Fibroblasts</i>	Adherent	90%	63%	
BAEC	<i>Bovine Aortic Endothelial Cells</i>	Adherent	80%	80%	Y

A6	<i>Xenopus Kidney Epithelial cells</i>	Adherent	90%	60%	Y
----	--	----------	-----	-----	---

## Primary Cells in Adherent Status: With Adherent-Cell Electrodes

Cell	Characteristics	V	TE	IMG
Human Skin Fibroblasts	Primary, Adherent	100%	50%	Y
HUVEC Human Umbilical Vein Endothelial cells	<i>Primary, Normal, Adherent</i>	95%	75%	Y
Mouse Cerebral Cortex Neurons (Embryonic day 15) (6 DIV)	Primary, Adherent	Excellent	Excellent	Y
Mouse Cerebral Cortex Neurons (Embryonic day 16) (2 DIV)	Primary, Adherent	70%	70%	Y
Mouse Cerebral Cortex Neurons (Embryonic day 16) (2 DIV)	Primary, Adherent	60%	65%	Y
Mouse Cerebral Cortex Neurons (Embryonic day 14) (5 DIV)	Primary, Adherent	95%	75%	Y
Mouse Hippocampal Neurons (Embryonic day 14) (4 DIV)	Primary, Adherent	60%	50%	
Mouse Hippocampal Neurons (Embryonic day 18) (2 DIV)	Primary, Adherent	85%	54%	Y
Mouse Neural Stem cells	Primary, Adherent	71%	50%	
Mouse Microglial cells (1 DIV after 1 week co-culturing astrocyte and microglial cells)	Primary, Adherent	Excellent	73%	Y
Mouse Glial cells (14 DIV)	Primary, Adherent	80%	50%	Y
Mouse Stromal cells (1-month cultured)	Primary Adherent	90%	50%	
Mouser Liver cells *siRNA Knock Down	Primary Adherent	Excellent	89%	
Rat Cerebral Cortex Neurons (Embryonic day 17) (2 DIV)	Primary, Adherent	70%	60%	Y
Rat Hippocampal Neurons (Embryonic day 16) (2 DIV)	Primary, Adherent	70%	70%	Y
Rat Hippocampal Neurons (Postnatal day 7) (11 DIV)	Primary, Adherent	100%	50%	Y
Rat Granulosa cells	Primary, Adherent	Excellent	41%	

## Cell Lines in Adherent Status: With Adherent-Cell Electrodes

Cell Line	Description, Species	Characteristics, etc.	V	TE	IMG
SH-SY5Y	<i>Human Neuroblastoma cells</i>	Adherent	90%	50%	Y
EPC	<i>Human Endothelial Progenitor Cells</i>	Adherent			Y
	<i>Human Mesenchymal Stem cells</i>	Bone Marrow derived	70%	65%	Y
THP-1	<i>Human Acute Monocytic Leukemia cells</i>	Suspension	90%	45%	
C2C12	<i>Mouse Myoblast cells</i>	Adherent	94%	>60%	Y
3T3-L1	<i>Mouse Embryonic Fibroblasts</i> <i>7 days after differentiation</i>	Adipocyte-like after differentiation	90%	70%	Y
MEF	<i>Mouse Embryonic Fibroblasts</i>	Adherent	60%	80%	
Neuro-2a	<i>Mouse Neuroblastoma cells</i>	Adherent	80%	90%	
C6	<i>Rat Glioma cells</i>	Fibroblast-like, GFAP Positive,	57%	55%	Y