


# in vitro zygote electroporation





Producing genetically modified mice by electroporation



- No microinjection
- Easier and faster
- Better efficiency

# in vitro zygote electroporation

## Electrode Types



### CUY501P1-1.5

- 1mm Gap
- 5 ul reagents
- 5 – 30 (or 50) zygotes
- Can reuse reagent

### CUY505P5

- 1mm Gap
- 50 ul reagents
- 20 – 150 zygotes
- Better results

# in vitro zygote electroporation

## Producing genetically modified mice by electroporation

Date	July 13, 2017						
Name	Dr. Tomoharu Yasuda						
Institution	Division of Immunology and Genomic Biology, Kyushu University						
Electroporator	<b>NEPA21 Typell</b>						
Electrode	5 mm gap electrode (CUY505P5)						
Electroporation buffer	1x OptiMEM						
Volume in electrode	45 ul (+/-α)						
Impedance	480-520 Ω						
Zygotes	C57BL/6J, IVF, Stored in fridge for 1 day						
Pre-electroporation wash	1x OptiMEM						
Post-electroporation wash/culture	KSOMaa						
Sample	Cas9 3NLS (IDT 1074181)	guide crRNA:tracrRNA (IDT)	Embryos electroporated	2-cell embryos (%)	Embryos transferred	Offspring (%)	Mutant/Offspring (%)
KO	100 ng/ul	36 ng/ul (1 uM)	88	54 (61%)	54	8 (15%)	8 (100%)
KO	100 ng/ul	36 ng/ul (1 uM)	88	52 (59%)	52	5 (10%)	5 (100%)
			<b>176</b>	<b>106 (60.2%)</b>	<b>106</b>	<b>13 (12.3%)</b>	<b>13 (100%)</b>

Data courtesy of Dr. Tomoharu Yasuda,  
Division of Immunology and Genomic Biology,  
Kyushu University



# in vitro zygote electroporation

## Efficient generation of knock-in rats using CRISPR/Cas9 Long KIs over 1500 bases

Rat embryos were collected from the superovulated females and cultured in a modified Krebs-Ringer bicarbonate medium (ARK resource, Kumamoto, Japan).

For EP, 400ng/ $\mu$ l Cas9 mRNA, or 100ng/ $\mu$ l Cas9 protein, 200 ng/ $\mu$ l gRNA, and 40ng/ $\mu$ l IssDNA were electroporated into rat embryos.

Rat embryos that developed to the two-cell stage after the introduction of RNA and IssDNA were transferred into the oviducts of female surrogates anesthetized with isoflurane (DS Pharma Animal Health Co., Ltd., Osaka, Japan).



**CONCLUSION** : Efficient production of KI rats by electroporation and Iss ODN is possible

Genes	KI lengths	Embryos electroporated	2cell embryos (%)	Embryos transferred	Pups delivered (%)	Mutant/offspring (%)	KI/offspring(%)
<i>GeneV</i> flox	674bp	160	77(48)	30	6(20)	6/6(100)	3/6(50)
<i>GeneG</i> flox	761bp	80	73(91)	60	27(45)	14/27(52)	2/14(7)
<i>GeneP</i> GFP	1421bp	180	131(72)	80	27(33)	27/27(100)	1/27(4)
<i>GeneM</i> 2A-iCre	1693bp	165	149(90)	60	11(18)	11/11(100)	1/11(9)

Data courtesy of Dr. Tomoji Mashimo,  
Genome Editing Research and Development Center,  
Graduate School of Medicine, Osaka University



# in vivo oviduct electroporation




## Genome-editing via Oviductal Nucleic Acids Delivery



Expose the oviduct of a pregnant mouse



Inject DNA into the ampulla of oviduct



Electroporation  
(transfer DNA into zygotes)

- **No embryo manipulation**
- **No pseudo-pregnant mice**
- **No vasectomized mice**
- **10-min treatment/mice**

Images courtesy of Dr. Shuji Takabayashi  
Laboratory for Animal Resources Development  
Hamamatsu University School Of Medicine



# in vivo oviduct electroporation

## Electrode Type



**CUY652P2.5X4 (old model No. CUY562)**

- **Concave platinum surface for maximum contact and minimum stress**

## in vivo oviduct electroporation

### Concentrations of Cas9 protein and efficiencies of GM-mouse production

CAS9 conc. (ng/ $\mu$ l)	No. of ICR mice	No. of Off- spring	Ratio (%) Offsprings /mice	No. of KI Offspring (%)	No. of Indel mutant (%)
50	11	107	9.7	2 ( 1.9)	4 ( 3.7)
100	7	47	6.7	20 (42.6)	23 (48.9)
500	6	45	7.5	28 (62.2)	32 (71.1)
1,000	7	47	6.7	28 (59.6)	36 (76.6)

Data courtesy of Dr. Shuji Takabayashi  
Laboratory for Animal Resources Development  
Hamamatsu University School Of Medicine

