

# Electrodes

## CUY560-5


### Parallel Fixed 2 Needle Straight Tip Electrode

10mm-Tip, 0.5mm diameter, 5mm gap

For animal tissues (Muscle, Skin of Mouse, Rat, Rabbit & Dog)

For DNA Vaccine

For Embryo, Testis, Bladder, Skin, Muscle and Organs



## CUY560-10, 15

### Parallel Fixed 2 Needle Straight Tip Electrode

10mm-Tip, 0.5mm diameter

CUY560-10: 10mm gap

CUY560-15: 15mm gap

For All animal tissues (Muscle, Skin of Rat, Rabbit, Dog & Pig)

For DNA Vaccine

